UREA NITROGEN

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GQ-BUN-60 1x60mL+ 1mL	GQ-BUN-100 1x100mL+ 1mL	GQ-BUN-125 1x125mL+ 1mL
Reagent used for the "In Vitro" quantitative determination of Urea Nitrogen in human serum.		

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PARAMETROS:

Temperature: 37°C
Wavelength: 340 nm
Assay type: Kinetic
Direction: decreasing

Ratio: Sample/Reagent 1:100
Sample volume: 0.003mL (3µL)
Reagent volume: 0.3
Tiempo de espera: 1 min
Lectura en min: 2 min

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Urea (mmol/L) = Urea (mg/dL) x 0.167

BIBLIOGRAPHIC REFERENCES

1. Zilva JF, Pannall PR. "UREA in Clinical Chemistry in Diagnosis and Treatment". Lloyd-Luke 1979; Chap XII; 286-8
2. Malloy HT, Evelyn KA. J Biol Chem 1937; 119:481-90
3. Jendrassik L, Grof B. Biochem Zeit 1938; 297; 81-9
4. Pearlman FC, Lee RT. Clin Chem 1974; 20:447-53
5. Young DS, Effects of Drugs on Clinical Laboratory Tests. Third Edition. AACC Press 1990
6. Tietz N, Textbook of Clinical Chemistry, 2nd Edition. Burtis CS and Ashwood ER (Eds) W.B. Saunders, 1995.
7. Watchel M et al, Creation and Verification of Reference Intervals. Laboratory Medicine, 1995;26:593-7

Manufactured by Bialex
717 Ohio Ave
Wichita Falls, TX 76301
United State of America

Contact us:

TollFree: 833-4182116

eFAX: +1-520-842-2001

Website: <http://www.bialexlab.com>

mail: info@bialexlab.com

To review the analysis protocol of your lot, visit:

Symbols use in the set

	Batch code	 YYYY-MM-DD	Useby/Expiration date
	Catalog Number		Manufacturer
	Authorized Representative		CE Mark
	In Vitro Diagnostic Medical Device		Temperature Limitation
	Consult Instructions for Use		Reagent 1/ Reagent 2

IFU-UREA BUN
Review 07-2019
Spanish version 1.0