



Composite Serum Control High – Level 3

Non Rutinary Biochemistry

Use: IVD
REF: 910726

Intended of Use

Composite Serum Control High Level was designed for use with different brands, since Bialex Lab uses other brands for QC, and compares the results with similar methods, however, it is highly recommended to use the same reagent—Bialex.

Preparation

Add 5mL of DI Water with pH 7.0 +/-0.2 to each vial using a precision volumetric pipette.

Let for 30 minutes or more, and from time to time stir gently.

After opened it is stable for 8 hours at 15 to 25°C and 3 days, at 2 to 8°C, and 04 weeks by -20°C.

Values:

Updates will be available at <http://platinum.bialexlab.com/files> Spanish language)
<http://platinum.bialexlab.com/> for English Version (U.S. Only)

Precaution

All specimens using in the Control should be handled in accordance with good laboratory practices using appropriate precautions as described in the CDC/NIH Manual, “Biosafety in Microbiological

and Biomedical Laboratories! 2nd edition. 1988 HHS Publication No (CDC)88-8395

DBA Bialex, certifies all individual human plasma donor units used to manufacture this lot were tested by FDA approved methods and found to ne nonreactive for HBsAg, HIV-1, HIV-2, and HCV

Control Serum for Non Rutinary Biochemistry (Level 3 – High values) LOT: 40413128UE Used
by: 2026.04.23

Anaylite	Value	Range High Level (C3)	Units
Bialex Cholinesterase	5278	3097-7459	U/L
Bialex CK	512	396-628	U/L
Bialex Urea UV	117 / 19.5	94-140 /15.6- 23.40	mg/dL or mmol/L
Bialex Chloride (REF 709141)	113 /113	101.7-124.3	mEq/L or mmol/L
Bialex Sodium (REF 71339)	162 /162	137.7-186.3	mEq/L or mmol/L
Bialex Ion	37.60	32-43.20	μmol/L
Bialex Total Bile Acid	44.3	42.09-46.51	μmol/L